

The development and application of a single cell biosensor for the detection of L-methionine and branched-chain amino acids

Nuriye Mustafi

Forschungszentrum Jülich GmbH
Institute of Bio- and Geosciences (IBG)
Biotechnology (IBG-1)

The development and application of a single cell biosensor for the detection of L-methionine and branched-chain amino acids

Nuriye Mustafi

Schriften des Forschungszentrums Jülich
Reihe Gesundheit / Health

Band / Volume 72


ISSN 1866-1785

ISBN 978-3-89336-956-0

Content

1	Summary	1
1.1.	Summary English	1
1.2.	Summary German	2
2	Introduction	4
2.1.	<i>Corynebacterium glutamicum</i> - A versatile microbial factory	4
2.2.	Biosynthesis of L-methionine in <i>C. glutamicum</i>	5
2.3.	Strategies for strain development	8
2.4.	Genetically-encoded biosensors	10
2.5.	Population heterogeneity in bacterial monocultures	12
2.6.	Aims of this work	14
3	Results	15
3.1.	Lrp of <i>Corynebacterium glutamicum</i> controls expression of the <i>brnFE</i> operon encoding the export system for L-methionine and branched-chain amino acids	18
3.2.	The development and application of a single-cell biosensor for the detection of L-methionine and branched-chain amino acids	30
3.3.	Destabilized eYFP variants for dynamic gene expression studies in <i>Corynebacterium glutamicum</i>	40
3.4.	L-methionine producing mutants of <i>Corynebacterium glutamicum</i> isolated by biosensor-based high throughput-screening	47
3.5.	Application of a genetically-encoded biosensor for the single cell analysis of <i>Corynebacterium glutamicum</i> L-valine production strains	56
3.6.	Development and application of genetically-encoded biosensors for strain development and single cell analysis of <i>Corynebacterium glutamicum</i>	74
4	Discussion	97
4.1.	Implementation of the Lrp-biosensor in FACS-based HT screenings	97
4.2.	Towards L-methionine overproduction	100

4.3.	Online monitoring of amino acid production	104
4.4.	Biosensor-based live cell imaging	104
4.5.	Limitations of the sensor system	107
5	References	109
6	Appendix	118
6.1.	Supplemental information (SI) - The development and application of a single-cell biosensor for the detection of L-methionine and branched-chain amino acids	118
6.2.	Supplemental information (SI) - L-methionine producing mutants of <i>Corynebacterium glutamicum</i> isolated by biosensor-based high throughput-screening	126
6.3.	Deletion of the aconitase gene in <i>Corynebacterium glutamicum</i> causes strong selection pressure for secondary mutations inactivating citrate synthase	128



Gesundheit / Health
Band / Volume 72
ISBN 978-3-89336-956-0

